



# FREZ ZE ZGRUBNYM OSTRZEM KRZYŻOWYM

Obróbka żywic i gipsu.



∅ mm / L

∅ mm / L

|                     |                     |                     |                     |
|---------------------|---------------------|---------------------|---------------------|
|                     |                     |                     |                     |
| 500 104 137 220 060 | 500 104 194 220 031 | 500 104 194 220 040 | 500 104 194 220 050 |
| 6,0 / 11,0          | 3,1 / 11,5          | 4,0 / 13,5          | 5,0 / 14,0          |

|                     |                     |                     |                     |
|---------------------|---------------------|---------------------|---------------------|
|                     |                     |                     |                     |
| 500 104 194 220 070 | 500 104 237 220 060 | 500 104 263 220 060 | 500 104 274 220 060 |
| 7,0 / 14,0          | 6,0 / 11,0          | 6,0 / 11,5          | 6,0 / 14,0          |





rodentica®

# FREZ Z OSTRZEM SPIRALNYM

Ostrze spiralne z poprzecznymi nacięciami  
Obróbka żywic.



500 104 194 176 040

500 104 225 176 060

500 104 274 176 060

500 104 277 176 040

∅ mm / L

40 / 13,5

6,0 / 12,5

6,0 / 14,0

4,0 / 8,0



[www.rodentica.eu](http://www.rodentica.eu)

Rodentica Polska Sp. z o.o.  
ul. Michejdy 12, 43-400 Cieszyn  
tel. : 388 412 428  
tel./fax: 033/8510717  
[www.facebook.com/rodentica](http://www.facebook.com/rodentica)